

# PLANT CELL WALL MONOCLONAL ANTIBODIES

>50 rat antibody hybridoma cell lines

Prof. Paul Knox

Centre for Plant Sciences, Faculty of Biological Sciences  
University of Leeds, Leeds LS2 9JT, UK

Email: [j.p.knox@leeds.ac.uk](mailto:j.p.knox@leeds.ac.uk)

[www.plantcellwalls.leeds.ac.uk](http://www.plantcellwalls.leeds.ac.uk)



## Notes on the selection and use of our MAbs

The LM- and other series of monoclonal antibodies (MAbs) developed in Leeds is a large resource of probes for the analysis of plant cell wall carbohydrates. If you require a specific antibody/epitope then that can be straightforward. If you wish to study and explore wider aspects of cell walls using a panel of antibodies that can detect a range of glycans then we have series of **recommended MAbs** for polysaccharide classes.

For pectic HG we recommend the combined use of **LM19** (unesterified HG) and **LM20** for high methylester HG. **JIM7** is also always additionally recommended for any analysis as it binds widely to pectin.

For detection of xyloglucan we recommend the high affinity **LM25** and for heteromannan **LM21**.

There is a range of probes for heteroxylan – the core two to use would be **LM11** that binds to 1,4-xylosyl residues and **LM28** to glucuronosyl residues. The often different binding patterns of **LM10** and **LM11** is now known to arise from the specific recognition by **LM10** of the non-reducing end (NRE) of xylans.

It can be hard to predict which AGP MAb to select as these glycan epitopes vary between tissues, organs and species. **JIM13** and **LM2** are a good place to start as they usually detect something in a section or an extract. Do not forget that a single AGP glycan epitope is unlikely to detect all AGPs in an organ. For starting with extensins we suggest the use of **LM1** and **JIM20**.

## Pectic polysaccharides

### Homogalacturonan (HG) / related

**LM19** partially Me-HG / no ester (34)

**LM20** partially Me-HG (34)

**JIM7** partially Me-HG (6,22,26)

**LM7** partially Me-HG / non-blockwise, (24,26)

**JIM5** partially Me-HG / no ester (1,6,22,26)

**LM18** partially Me-HG / no ester (34)

**PAM1** blockwise de-esterified HG (17,22,23,30)

**LM8** xylogalacturonan (27)

### Rhamnogalacturonan-I

#### Galactan +

**LM5 NRE** (1→4)-β-D-galactan (15,18-20,42)

**LM26** branched (1,6-Gal) (1→4)-β-D-galactan (45)

#### Arabinan +

**LM6<sup>+</sup>** (1→5)-α-L-arabinan (16,18-20,23,31,34)

**LM6-M** (1→5)-α-L-arabinan (46)

**LM13** linearised (1→5)-α-L-arabinan (33,35)

#### Other

**LM16** processed arabinan/put. galactan stub (35)

**LM9** feruloylated (1→4)-β-D-galactan (28)

**LM12<sup>‡</sup>** ferulic acid, feruloylated pectin (38)

<sup>+</sup>May also bind to AGPs

<sup>‡</sup>Can also bind to feruloylated heteroxylan

## Non-cellulosic, non-pectic polysaccharides

### Xyloglucan

- LM15** XXXG motif of xyloglucan (32,44)
- LM24** galactosylated xyloglucan (38)
- LM25** XXXG/galactosylated xyloglucan (38)

### Heteromannan

- LM21** heteromannan (36)
- LM22** heteromannan (36)

### Heteroxylan

- LM10 NRE** (1→4)-β-D-xylan (29,44)
- LM11** (1→4)-β-D-xylan / arabinoxylan (29)
- LM28** glucuronoxylan (40)
- LM12<sup>‡</sup>** ferulic acid, feruloylated xylan (38)
- LM27** unknown epitope assoc. grass xylan (40)

<sup>‡</sup>Can also bind to feruloylated pectin

## Plant cell wall proteoglycans/glycoproteins

### Arabinogalactan-protein (AGP) glycan

- LM2** β-linked-GlcA in AGP glycan (13,14)
- LM14** GlcA in AGP glycan (33,38)
- LM30** AGP glycan (43)
- JIM4** AGP glycan (3,5,13)
- JIM13** AGP glycan (7,13)
- JIM14** AGP glycan (7,13,44)
- JIM15** AGP glycan (7,13)
- JIM16** AGP glycan (7,13,44)
- MAC207** AGP glycan (2,3,13)

### Extensin

- LM1** extensin (11)
- JIM11** extensin (8)
- JIM12** extensin (8)
- JIM19** extensin (8,9,10)
- JIM20** extensin (8,9)

**NRE** = epitope at non-reducing end of glycan

## Other cell wall related MAbs

- LM23** non-acetylated xylosyl in xylogalacturonan, xylan, fucoidan preps (37,38, 39)
- LM4** pea amine oxidase, cell walls (21,25)
- JIM18** glyco-phospholipid, membranes (9,12)
- JIM1** β-linked-galactosyl, plasma membrane (4)

## Brown algal cell wall polysaccharides

### Fucoidan

- BAM1** un-sulfated epitope present in sulfated fucan/fucoidan preparations (39)
- BAM2** sulfated epitope present in sulfated fucan/fucoidan preparations (39)
- BAM3** possibly sulfated epitope present in sulfated fucan/fucoidan preparations (39)
- BAM4** sulfated epitope present in sulfated fucan/fucoidan preparations (39)

### Alginate

- BAM6** mannuronate-rich epitope (41)
- BAM7** mannuronate-guluronate (41)
- BAM8** mannuronate-guluronate (41)
- BAM9** mannuronate-guluronate (41)
- BAM10** mannuronate-guluronate epitope resistant to alginate lyase (41)
- BAM11** ~7 guluronate residues (41)

If you have specific questions about our MAb specificities or use please enquire at [j.p.knox@leeds.ac.uk](mailto:j.p.knox@leeds.ac.uk)

## References

1. VandenBosch *et al.* (1989) **EMBO Journal** 8, 335-342
2. Pennell. *et al.* (1989) **J. Cell Biology** 108, 1967-1977
3. Knox *et al.* (1989) **Development** 106, 47-56
4. Knox & Roberts (1989) **Protoplasma** 152, 123-129
5. Stacey *et al.* (1990) **Planta** 180, 285-292
6. Knox *et al.* (1990) **Planta** 181, 512-521
7. Knox, *et al.* (1991) **Plant Journal** 1, 317-326
8. Smallwood *et al.* (1994) **Plant Journal** 5, 237-246
9. Knox *et al.* (1995) **Planta** 196, 266-270
10. Wang *et al.* (1995) **Planta** 196, 271-276
11. Smallwood *et al.* (1995) **Planta** 196, 510-522
12. Perotto *et al.* (1995) **Mol. Plant-Microbe Int.** 8, 560-568
13. Yates *et al.* (1996) **Glycobiology** 6, 131-139
14. Smallwood *et al.* (1996) **Planta** 198, 452-459
15. Jones *et al.* (1997) **Plant Physiology** 113, 1405-1412
16. Willats *et al.* (1998) **Carbohydr. Research** 308, 149-152
17. Willats *et al.* (1999) **Plant Journal** 18, 57-65
18. Willats *et al.* (1999) **Plant Journal** 20, 619-628
19. McCartney *et al.* (2000) **Plant Journal** 22, 105-113
20. Orfila & Knox (2000) **Plant Physiology** 122, 775-781
21. Wisniewski *et al.* (2000) **Mol. Plant-Microbe Int.** 13, 413
22. Willats *et al.* (2000) **Carbohydr. Res.** 327, 309-320
23. Orfila *et al.* (2001) **Plant Physiology** 126, 210-221
24. Willats *et al.* (2001) **J. Biol. Chem.** 276, 19404-19413
25. Laurenzi *et al.* (2001) **Planta** 214, 37-45
26. Clausen *et al.* (2003) **Carbohydr. Res.** 338, 1797-1800
27. Willats *et al.* (2004) **Planta** 218, 673-681
28. Clausen *et al.* (2004) **Planta** 219, 1036-1041
29. McCartney *et al.* (2005) **J. Histochem Cytochem** 53, 543
30. Manfield *et al.* (2005) **Plant Science** 169, 1090-1095
31. Lee *et al.* (2005) **Plant Cell** 17, 3051-3065
32. Marcus *et al.* (2008) **BMC Plant Biology** 8, 60
33. Moller *et al.* (2008) **Glycoconjugate J.** 25, 37-48
34. Verherbruggen *et al.* (2009) **Carbohydr. Res.** 344, 1858
35. Verherbruggen *et al.* (2009) **Plant Journal** 59, 413-425
36. Marcus *et al.* (2010) **Plant Journal** 64, 191-203
37. Manabe *et al.* (2011) **Plant Physiology** 155, 1068-1078
38. Pedersen *et al.* (2012) **J. Biol. Chem.** 287, 39429-39438
39. Torode *et al.* (2015) **PlosONE** 10, e0118366
40. Cornuault *et al.* (2015) **Planta** 242, 1321-1334
41. Torode *et al.* (2016) **J. Exp Botany** 67, 6089-6100
42. Andersen *et al.* (2016) **Carbohydrate Research** 436, 36-40
43. Wilkinson *et al.* (2017) **J. Cereal Science** 74, 155-164
44. Ruprecht *et al.* (2017) **Plant Physiology** 175, 1094-1104
45. Torode *et al.* (2018) **Plant Physiology** 176, 1547-1558
46. Cornuault *et al.* (2018) **BioRxiv** /10.1101/161604

## MAB selections available through

### Kerafast

<https://www.kerafast.com/cat/799/paul-knox-phd>

### Megazyme

<https://www.megazyme.com/shop-all-products/antibodies>

### Ximbio

<https://ximbio.com/search?q=paul+knox>

## Brown algal BAM-series of MABs available from SeaProbes

<http://www.sb-roscoff.fr/en/seaprobes>

Autumn 2021